

## Product datasheet for **MP203384**

### Dcun1d1 Mouse qPCR Primer Pair (NM\_033623)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                                   |
|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                                 |
| Gene ID:          | 114893                                                                                                                                                                                                                                                                                                                                            |
| Forward Sequence: | GTTGGACAGGAAGAAGTTAGAGC                                                                                                                                                                                                                                                                                                                           |
| Reverse Sequence: | TCAACACACTGATGCTGGCTGG                                                                                                                                                                                                                                                                                                                            |
| ACCN:             | <u>BC020161</u> , <u>BC031666</u> , <u>BC037431</u> , <u>NM_033623</u> , <u>NM_033623.1</u> , <u>NM_033623.2</u> , <u>NM_033623.3</u> ,<br><u>NM_033623.4</u> , <u>NM_033623.5</u> , <u>NM_033623.6</u> , <u>BC020161.1</u> , <u>BB751776</u> , <u>BY392516</u>                                                                                   |
| UniProt ID:       | <u>Q9QZ73</u>                                                                                                                                                                                                                                                                                                                                     |
| Synonyms:         | pTes3; Rp42; SCCRO; Tes3                                                                                                                                                                                                                                                                                                                          |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions).<br>Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                    |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec. |
| Storage:          | Store at -20 °C.                                                                                                                                                                                                                                                                                                                                  |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                      |

