

## Product datasheet for **MP203388**

### Ddx17 Mouse qPCR Primer Pair (NM\_199080)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67040
Forward Sequence:	GACAAAGAGGCGCTGTGATGAC
Reverse Sequence:	AGGAGCCTTTCCAGATCGGAAC
ACCN:	<a href="#">BC062910</a> , <a href="#">BC096036</a> , <a href="#">NM_199080</a> , <a href="#">NM_199080.1</a> , <a href="#">NM_199080.2</a> , <a href="#">BC012249</a> , <a href="#">BC027758</a> , <a href="#">BC031802</a> , <a href="#">BC038378</a> , <a href="#">BI655302</a>
UniProt ID:	<a href="#">Q501J6</a>
Synonyms:	2610007K22Rik; A430025E01Rik; AI047725; C80929; Gm926; p7; p72
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)