

Product datasheet for **MP203389**

Ddx18 Mouse qPCR Primer Pair (NM_025860)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66942
Forward Sequence:	CTCGCTTTCCTCATCCCTGTGA
Reverse Sequence:	GTGTGAACGTGGTGCATCA
ACCN:	BC028246 , BC103776 , NM_025860 , NM_025860.1 , NM_025860.2 , NM_025860.3 , BC017160 , BC141230 , BC141231 , BY084565 , BY381480
UniProt ID:	Q8K363
Synonyms:	2310005B10Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)