

Product datasheet for **MP203397**

Ddx27 Mouse qPCR Primer Pair (NM_153065)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	228889
Forward Sequence:	AGCACCAGACATCCTCATTGCC
Reverse Sequence:	GCTCCTCGAAATACTCATCCAGC
ACCN:	<u>BC011321</u> , <u>NM_153065</u> , <u>NM_153065.1</u> , <u>NM_153065.2</u> , <u>NM_153065.3</u> , <u>BC024730</u> , <u>BC024981</u> , <u>BC026381</u> , <u>BC033398</u>
UniProt ID:	<u>Q92IN6</u>
Synonyms:	C86129
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

