

Product datasheet for **MP203403**

Ddx42 Mouse qPCR Primer Pair (NM_028074)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72047
Forward Sequence:	CAGGTTTCGGTCCATAGCAAGTC
Reverse Sequence:	ATGTCTCCCTGCACCACTCGAA
ACCN:	BC043036 , NM_028074 , NM_028074.1 , NM_028074.2 , NM_028074.3 , BC019772 , BC049382 , BC055482
UniProt ID:	Q810A7
Synonyms:	1810047H21Rik; AW319508; AW556242; B430002H05Rik; RHELP; RNAHP; SF3b125
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)