

Product datasheet for **MP203404**

Ddx46 Mouse qPCR Primer Pair (NM_145975)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	212880
Forward Sequence:	CCTTGCATGTCTCTTCACGGAG
Reverse Sequence:	ACGAGCTGCAACAGATGTGGCA
ACCN:	BC026492 , BC092240 , NM_145975 , NM_145975.1 , NM_145975.2 , NM_145975.3 , BC026492.1 , BC019403 , BC028552
UniProt ID:	Q569Z5
Synonyms:	2200005K02Rik; 8430438J23Rik; AI325430; AI957095; mKIAA0801
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)