

Product datasheet for **MP203429**

Dhx30 Mouse qPCR Primer Pair (NM_133347)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72831
Forward Sequence:	CCGTCAAGGAACACTACCTGGA
Reverse Sequence:	CAGAACCAGGTCAGTCACAAGG
ACCN:	<u>BC016202</u> , <u>NM_133347</u> , <u>NM_133347.1</u> , <u>NM_133347.2</u> , <u>BC016202.1</u> , <u>BC004082</u> , <u>BC058331</u> , <u>BC062885</u> , <u>BC069884</u> , <u>BC084729</u> , <u>BY540912</u>
UniProt ID:	<u>Q99PU8</u>
Synonyms:	2810477H02Rik; C130058C04Rik; Ddx30; HELG; Ret-CoR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

