

Product datasheet for **MP203431**

Dhx33 Mouse qPCR Primer Pair (NM_178367)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	216877
Forward Sequence:	TGGCAGCATTCCCTTTAGAGCC
Reverse Sequence:	GACACTGTCCACAGATAGCAGG
ACCN:	<u>BC052172</u> , <u>NM_178367</u> , <u>NM_178367.2</u> , <u>NM_178367.3</u> , <u>NM_178367.4</u> , <u>BC011309</u> , <u>BC027842</u> , <u>BQ287292</u>
UniProt ID:	<u>Q80VY9</u>
Synonyms:	3110057P17Rik; 9430096J02Rik; Ddx33
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

