

Product datasheet for **MP203439**

Dhx9 Mouse qPCR Primer Pair (NM_007842)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13211
Forward Sequence:	GGCAATCGAACCTCCACCTTTG
Reverse Sequence:	GGAGTTTAGCCAGGATTCGTCC
ACCN:	<u>NM_007842</u> , <u>NM_007842.1</u> , <u>NM_007842.2</u> , <u>BC033911</u> , <u>BC064082</u> , <u>BC089159</u> , <u>BC169284</u> , <u>BC169285</u> , <u>BC169286</u>
UniProt ID:	<u>O70133</u>
Synonyms:	AI326842; Ddx9; HEL-5; mHEL-5; NDHII; RHA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

