

Product datasheet for **MP203442**

Dapk1 Mouse qPCR Primer Pair (NM_029653)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69635
Forward Sequence:	AGATGTGGTCCGCTACCTCTGT
Reverse Sequence:	ATGCTCGTGCTGTTCTGCCTTG
ACCN:	<u>BC057317</u> , <u>BC060161</u> , <u>NM_029653</u> , <u>NM_029653.1</u> , <u>NM_029653.2</u> , <u>NM_029653.3</u> , <u>BC021490</u> , <u>BC026671</u> , <u>BC037483</u> , <u>BQ175948</u> , <u>BQ745486</u> , <u>BY093371</u>
UniProt ID:	<u>Q80YE7</u>
Synonyms:	D13Ucla1; DAP-Kinase
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)