

Product datasheet for **MP203453**

Dcn Mouse qPCR Primer Pair (NM_007833)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13179
Forward Sequence:	ACTCTCCAGGAAGCTTCGTGTCC
Reverse Sequence:	AGTCCCTGGAAGGCTCCGTTTT
ACCN:	BC132521 , NM_007833 , NM_007833.1 , NM_007833.2 , NM_007833.3 , NM_007833.4 , NM_007833.5 , NM_007833.6 , BB575065 , BC060126 , BC138564 , BG088078 , BM114383
UniProt ID:	P28654
Synonyms:	DC; DSPG2; PG40; PGII; PGS2; SL; SLRR1B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)