

Product datasheet for **MP203459**

Dock6 Mouse qPCR Primer Pair (NM_177030)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	319899
Forward Sequence:	TGTTGACTCCCGCTACGCAGAA
Reverse Sequence:	ATGGAAGCCAGTCTTGACCGCT
ACCN:	<u>BC057368</u> , <u>NM_177030</u> , <u>NM_177030.1</u> , <u>NM_177030.2</u> , <u>NM_177030.3</u> , <u>BC021414</u> , <u>BC037714</u> , <u>BC043042</u> , <u>BU055715</u> , <u>BU701914</u> , <u>BU705495</u>
UniProt ID:	<u>Q8VDR9</u>
Synonyms:	2410095B20Rik; 4931431C02Rik; C330023D02Rik; mKIAA1395
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

