

Product datasheet for **MP203460**

Dock7 Mouse qPCR Primer Pair (NM_026082)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67299
Forward Sequence:	GAGAGACTCAACACAAGCCTTGC
Reverse Sequence:	CCAGGACACTTGGATTTGGCAG
ACCN:	<u>BC006868</u> , <u>BC086672</u> , <u>NM_026082</u> , <u>NM_026082.1</u> , <u>NM_026082.2</u> , <u>NM_026082.3</u> , <u>NM_026082.4</u> , <u>NM_026082.5</u> , <u>BC018599</u> , <u>BC024823</u> , <u>BC024917</u> , <u>BC055032</u> , <u>BC094420</u> , <u>BF141712</u> , <u>BG067751</u> , <u>BM964053</u> , <u>BQ938873</u> , <u>BY710986</u>
Synonyms:	3110056M06Rik; Gm430
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

