

Product datasheet for **MP203461**

Dock8 Mouse qPCR Primer Pair (NM_028785)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	76088
Forward Sequence:	TGCTCTGGCAATAGTGATGCGC
Reverse Sequence:	CTCACTACTCCAGTGCTGACCA
ACCN:	<u>NM_028785</u> , <u>NM_028785.3</u> , <u>BC141358</u> , <u>BC029018</u> , <u>BC030316</u> , <u>BC043470</u> , <u>BC055295</u>
UniProt ID:	<u>Q8C147</u>
Synonyms:	I200017A24Rik; 5830472H07Rik; A130095G14Rik; A1461977
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

