

Product datasheet for **MP203515**

Defa-rs2 Mouse qPCR Primer Pair (NM_007847)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13222
Forward Sequence:	CTAAAACTGAGGAGCAGCCAGG
Reverse Sequence:	ATGACGGGCATCTAGGGCATTG
ACCN:	<u>BC132646</u> , <u>BC132648</u> , <u>NM_007847</u> , <u>NM_007847.1</u>
UniProt ID:	<u>P17534</u>
Synonyms:	CRS4C; CRS4C-1; CRS4C1d; CRS4C1e; CRS4C1f; CRS4C1g; CRS4C1h; CRS4C1i; CRS4C1j; Defcr-rs2; Defcr-rs3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)