

Product datasheet for **MP203517**

Defa29 Mouse qPCR Primer Pair (NM_007844)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13218
Forward Sequence:	GTCTCCTTTGGAGGCACAGAAG
Reverse Sequence:	CACAATGGCTGCTTTGGGCACT
ACCN:	<u>BC117042</u> , <u>BC117044</u> , <u>NM_007844</u> , <u>NM_007844.1</u> , <u>NM_007844.2</u> , <u>BC145556</u> , <u>BC171966</u>
Synonyms:	2010300L12Rik; 2010319H24Rik; CRS1C; CRS1C-2; CRS1C-4; CRS1C-5; Defa-rs1; Defcr-rs1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

