

Product datasheet for **MP203541**

Dlc1 Mouse qPCR Primer Pair (NM_015802)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50768
Forward Sequence:	ATGACCGAAGGACACCCAGTGA
Reverse Sequence:	AAGCTGTGCCACCTCAGTCTGT
ACCN:	NM_015802 , NM_015802.1 , NM_015802.2 , NM_015802.3 , BC150802 , BC028562 , BC048730 , BC056230 , BC096052 , BE199599
UniProt ID:	Q9R0Z9
Synonyms:	A730069N07Rik; Arhgap7; dlc-1; HP; STARD12
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)