

Product datasheet for **MP203568**

Dennd4c Mouse qPCR Primer Pair (NM_001081014)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	329877
Forward Sequence:	CTGAACCTCAGAGTCCTACGGA
Reverse Sequence:	TGCATCCTCGACTGGAGCTTGA
ACCN:	<u>NM_001081014</u> , <u>NM_001081014.1</u> , <u>BC040249</u> , <u>BC040798</u> , <u>BC146602</u>
UniProt ID:	<u>A6H8H2</u>
Synonyms:	I700065A05Rik; AA420392; RP24-468M3.1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

