

Product datasheet for **MP203589**

Depdc1b Mouse qPCR Primer Pair (NM_178683)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	218581
Forward Sequence:	GAGTGGTCATCCTGGACGACAA
Reverse Sequence:	CATAGGCTGCTTCCAGTCTGTAC
ACCN:	<u>NM_178683</u> , <u>NM_178683.1</u> , <u>NM_178683.2</u> , <u>NM_178683.3</u> , <u>NM_178683.4</u> , <u>BC152330</u> , <u>BC048246</u> , <u>BC152331</u>
UniProt ID:	<u>Q8BH88</u>
Synonyms:	9830132O03; AW260467; XTP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

