

Product datasheet for **MP203593**

Dcakd Mouse qPCR Primer Pair (NM_026551)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68087
Forward Sequence:	CCGATACGTGATTCTGGACATCC
Reverse Sequence:	CAGGTTGTTCCGCTTCATCAGC
ACCN:	<u>BC005607</u> , <u>BC045193</u> , <u>NM_026551</u> , <u>NM_026551.2</u> , <u>NM_026551.3</u> , <u>BC023208</u> , <u>BC031855</u> , <u>BC042685</u> , <u>BC047340</u> , <u>BC058591</u>
UniProt ID:	<u>Q8BHC4</u>
Synonyms:	3010024O21Rik; 6720485C15Rik; A1849483
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

