

## Product datasheet for **MP203654**

### Dicer1 Mouse qPCR Primer Pair (NM\_148948)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	192119
Forward Sequence:	AGCAGTGCTGAGAAGAGGAAGG
Reverse Sequence:	CCGCTTTTCTCCACAGTGATGC
ACCN:	<u><a href="#">NM_148948</a></u> , <u><a href="#">NM_148948.1</a></u> , <u><a href="#">NM_148948.2</a></u> , <u><a href="#">BC019805</a></u> , <u><a href="#">BC023863</a></u> , <u><a href="#">BC029665</a></u> , <u><a href="#">BC061198</a></u> , <u><a href="#">BC100374</a></u> , <u><a href="#">BC172621</a></u>
Synonyms:	1110006F08Rik; D12Ert7e; mKIAA0928
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)