

## Product datasheet for **MP203700**

### Dpp6 Mouse qPCR Primer Pair (NM\_010075)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13483
Forward Sequence:	CAATCCTGACCCTCTGTGATGC
Reverse Sequence:	GCCATCTTTGGAGAACACAGGC
ACCN:	<u><a href="#">BC048383</a></u> , <u><a href="#">BC085154</a></u> , <u><a href="#">NM_010075</a></u> , <u><a href="#">NM_010075.1</a></u> , <u><a href="#">NM_010075.2</a></u> , <u><a href="#">BY127642</a></u> , <u><a href="#">BY280444</a></u>
UniProt ID:	<u><a href="#">Q9Z218</a></u>
Synonyms:	B930011P16Rik; D5Buc3; D5Buc4; D5Buc5; Dpp-6; DPP VI; Gm1377; In(5)6H-p; Peplb; Rw
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)