

Product datasheet for **MP203759**

Dna2 Mouse qPCR Primer Pair (NM_177372)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	327762
Forward Sequence:	GGAACCTCAGTTCATCGCCTAC
Reverse Sequence:	TGGCATCCCAACAATGAGCGTG
ACCN:	<u>BC115716</u> , <u>NM_177372</u> , <u>NM_177372.1</u> , <u>NM_177372.2</u> , <u>NM_177372.3</u> , <u>BC025182</u> , <u>BC030382</u> , <u>BM202919</u> , <u>BY768981</u>
UniProt ID:	<u>Q6ZQJ5</u>
Synonyms:	Dna2l; E130315B21Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)