

Product datasheet for **MP203766**

Mrpl42 Mouse qPCR Primer Pair (NM_026065)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67270
Forward Sequence:	ACAGCCGACGGCAGGACAATAG
Reverse Sequence:	TGGCTTTCAGGATCTGCTCGTG
ACCN:	<u>BC024337, NM_026065, NM_026065.1, NM_026065.2, NM_026065.3, NM_026065.4</u>
UniProt ID:	<u>Q9CPV3</u>
Synonyms:	2700009F22Rik; 2900055D03Rik; D10ErtD322e; HSPC204; L31mt; MRP-L31; PTD007; Rplm131
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)