

Product datasheet for **MP203789**

D16Ertd472e Mouse qPCR Primer Pair (NM_025967)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67102
Forward Sequence:	CCAGAATCCACATTACAGCCGTG
Reverse Sequence:	GAGAGTCCTTCTTGCACCTG
ACCN:	<u>BC019957, NM_025967, NM_025967.1, NM_025967.2, NM_025967.3, NM_025967.4, BG078224</u>
UniProt ID:	<u>Q9D7G4</u>
Synonyms:	1700010I10Rik; 2310009O17Rik; E330003K22Rik; Eurl
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

