

Product datasheet for **MP203797**

Pum3 Mouse qPCR Primer Pair (NM_177474)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	52874
Forward Sequence:	GGTCATGTGAGGAAGATGCTGC
Reverse Sequence:	TTGTCCAGCGTTGGGTGATCTG
ACCN:	<u>NM_177474</u> , <u>NM_177474.1</u> , <u>NM_177474.2</u> , <u>NM_177474.3</u> , <u>NM_177474.4</u> , <u>NM_177474.5</u> , <u>BC016186</u> , <u>BC041111</u> , <u>BC052686</u> , <u>BC145533</u> , <u>BY251590</u> , <u>BY289042</u>
UniProt ID:	<u>Q8BKS9</u>
Synonyms:	1110069H02Rik; AA675048; D19Bwgl357e
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

