

Product datasheet for **MP203806**

Plppr4 Mouse qPCR Primer Pair (NM_177664)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	229791
Forward Sequence:	CTTGGCTTTTGCTGGACCTGCA
Reverse Sequence:	GGCGTTGATGTTAGGCTCCAATC
ACCN:	<u>NM_177664</u> , <u>NM_177664.1</u> , <u>NM_177664.2</u> , <u>NM_177664.3</u> , <u>NM_177664.4</u> , <u>NM_177664.5</u> , <u>BC137702</u> , <u>BC024711</u> , <u>BC082323</u> , <u>BC137701</u>
UniProt ID:	<u>Q7TME0</u>
Synonyms:	A330086D10; D3Bwg0562e; Lppr4; mKIAA0455; PRG-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

