

Product datasheet for **MP203818**

Selenoi Mouse qPCR Primer Pair (NM_027652)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	28042
Forward Sequence:	GCTTTGGATACCAACCCACTCTC
Reverse Sequence:	GGTCGAAGTATGTCAGGAGTAGG
ACCN:	<u>BC106097</u> , <u>BC115783</u> , <u>BC117551</u> , <u>NM_027652</u> , <u>NM_027652.1</u> , <u>NM_027652.2</u> , <u>NM_027652.3</u> , <u>BC051183</u> , <u>BC055810</u> , <u>BC078438</u> , <u>BY136287</u>
UniProt ID:	<u>Q80TA1</u>
Synonyms:	4933402G07Rik; AI448296; AI452230; C79563; D5Wsu178; D5Wsu178e; Ept; Eptl; mKIAA1724; SELI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

