

Product datasheet for **MP203831**

Nrep Mouse qPCR Primer Pair (NM_053078)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27528
Forward Sequence:	GGTGTCGGTACTTTGTTTCCTGG
Reverse Sequence:	CTCACACTCTTGGTAGCATCCAC
ACCN:	BC054762 , NM_053078 , NM_053078.1 , NM_053078.2 , NM_053078.3 , NM_053078.4 , BQ715224 , BU057200 , BY150793
UniProt ID:	Q07475
Synonyms:	AI325076; D0H4S114; Harp; P311; PTZ17; SEZ17
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)