

Product datasheet for **MP203880**

Dok7 Mouse qPCR Primer Pair (NM_172708)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	231134
Forward Sequence:	TCAGCCTCAGAAGAGCGTGTG
Reverse Sequence:	GCCTCAGAAGAGGAAGTGGATAG
ACCN:	<u>BC089590</u> , <u>NM_172708</u> , <u>NM_172708.1</u> , <u>NM_172708.2</u> , <u>NM_172708.3</u> , <u>BC150489</u> , <u>NM_172708.4</u>
Synonyms:	A930013K19Rik; AW049091; Dok-7; EF-12; Oit5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)