

Product datasheet for **MP203925**

Dph2 Mouse qPCR Primer Pair (NM_026344)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67728
Forward Sequence:	CCACCTCCTTGGAAGTCATCAG
Reverse Sequence:	GACGGCTTCTTCACTGGTGTC
ACCN:	<u>BC031124</u> , <u>NM_026344</u> , <u>NM_026344.1</u> , <u>NM_026344.2</u> , <u>NM_026344.3</u> , <u>BC014746</u>
UniProt ID:	<u>Q9CR25</u>
Synonyms:	9130020C19Rik; AI467389; Dph2l2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)