

Product datasheet for **MP203977**

Dym Mouse qPCR Primer Pair (NM_027727)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69190
Forward Sequence:	CCTTCTGGAAGAGTTGCTCTGC
Reverse Sequence:	ACCTCCTTGTGGAAGAGTTGGC
ACCN:	<u>BC018220</u> , <u>BC038276</u> , <u>NM_027727</u> , <u>NM_027727.1</u> , <u>NM_027727.2</u> , <u>BC018220.1</u> , <u>BC092283</u> , <u>NM_027727.3</u>
UniProt ID:	<u>Q8CHY3</u>
Synonyms:	1810041M12Rik; 4933427L07Rik; C030019K18Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

