

Product datasheet for **MP203987**

Dnm3 Mouse qPCR Primer Pair (NM_001038619)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	103967
Forward Sequence:	GTTGCTAACCACATCCGTGAGAG
Reverse Sequence:	CTTGTGAACCTGACTGCTCCTC
ACCN:	<u>NM_001038619</u> , <u>NM_001038619.1</u> , <u>BC141143</u> , <u>BC004050</u> , <u>BC024584</u> , <u>BC141144</u> , <u>NM_001038619.2</u>
UniProt ID:	<u>Q8BZ98</u>
Synonyms:	9630020E24Rik; AW061159; B230343F03Rik; C530045C17; mKIAA0820
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

