

Product datasheet for **MP203996**

DynII2 Mouse qPCR Primer Pair (NM_026556)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68097
Forward Sequence:	GCCGTTGACTGTGCCACACAG
Reverse Sequence:	GTCTCGTGTGTGACATAGCTGC
ACCN:	<u>BC011289</u> , <u>BC040822</u> , <u>NM_026556</u> , <u>NM_026556.1</u> , <u>NM_026556.2</u> , <u>NM_026556.3</u> , <u>NM_026556.4</u> , <u>BC049621</u> , <u>BC058797</u> , <u>BY249698</u> , <u>BY310878</u>
UniProt ID:	<u>Q9D0M5</u>
Synonyms:	1700064A15Rik; 6720463E02Rik; C87222; Dlc2; DLC8; DLC8b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)