

Product datasheet for **MP204008**

Dnal1 Mouse qPCR Primer Pair (NM_028821)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	105000
Forward Sequence:	ATTGAGAAGTTGAAAGGAATCCATG
Reverse Sequence:	GGATTACCGACGAACACCAGGT
ACCN:	BC125392 , BC125394 , NM_028821 , NM_028821.1 , NM_028821.2 , NM_028821.3 , BC070454 , BY126068
UniProt ID:	Q05A62
Synonyms:	1700010H15Rik; AW121714; Dnalc1; E330027P08Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)