

Product datasheet for **MP204019**

Dtna Mouse qPCR Primer Pair (NM_207650)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13527
Forward Sequence:	CGGCTTGATGAAGAACACAGGC
Reverse Sequence:	CGATGGTGAAGGAGATGTCAGG
ACCN:	<u>BC040364</u> , <u>BC050765</u> , <u>NM_207650</u> , <u>NM_207650.1</u> , <u>NM_207650.2</u> , <u>NM_207650.3</u> , <u>NM_207650.4</u> , <u>BB337232</u> , <u>BU700818</u> , <u>BY122991</u> , <u>BY127150</u> , <u>BY130082</u>
UniProt ID:	<u>Q9D2N4</u>
Synonyms:	2210407P21Rik; a-DB-1; adbn; Dtn; DTN-A; Gm19389
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)