

Product datasheet for **MP204030**

Ets1 Mouse qPCR Primer Pair (NM_011808)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23871
Forward Sequence:	CCAGAATCCTGTTACACCTCGG
Reverse Sequence:	CAGCGTCTGATAGGACTCTGTG
ACCN:	<u>BC010588</u> , <u>BC114351</u> , <u>NM_011808</u> , <u>NM_011808.1</u> , <u>NM_011808.2</u> , <u>BB621050</u> , <u>BC013089</u>
UniProt ID:	<u>P27577</u>
Synonyms:	AI196000; AI448617; D230050P06; Ets-1; Tpl1; vs
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)