

Product datasheet for **MP204043**

Elf2 Mouse qPCR Primer Pair (NM_023502)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69257
Forward Sequence:	CACAGAGACTGTAGTGGAGGTG
Reverse Sequence:	GGAGACCCATTGGAACTGGTG
ACCN:	<u>BC092057</u> , <u>NM_023502</u> , <u>NM_023502.1</u> , <u>NM_023502.2</u> , <u>BC005619</u> , <u>BC027739</u> , <u>BC051937</u> , <u>BC069901</u>
UniProt ID:	<u>Q9JHC9</u>
Synonyms:	2610036A20Rik; A230104O07Rik; AW111824; BB183398; EU32; NERF; NERF-1A; NERF-1B; NERF-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)