

## Product datasheet for **MP204078**

### Entpd4 Mouse qPCR Primer Pair (NM\_026174)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67464
Forward Sequence:	TGGGTTTGGTGGTAATGCTGCC
Reverse Sequence:	TCCAGTAGTGGAGCATCAGGAG
ACCN:	<a href="#">BC043134</a> , <a href="#">NM_026174</a> , <a href="#">NM_026174.1</a> , <a href="#">NM_026174.2</a> , <a href="#">NM_026174.3</a> , <a href="#">BC006924</a> , <a href="#">BC026420</a> , <a href="#">BC052090</a> , <a href="#">BY276189</a>
UniProt ID:	<a href="#">Q9DBT4</a>
Synonyms:	1200014F22Rik; 4930513O14Rik; A430102M06Rik; AI414893; Lysal1; mKIAA0392
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)