

Product datasheet for **MP204107**

Egfl7 Mouse qPCR Primer Pair (NM_178444)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	353156
Forward Sequence:	CTGAGCATGTCTACAGACCCAG
Reverse Sequence:	TCTGTGTCCGTCGCAAGTGGTG
ACCN:	<u>BC024610</u> , <u>NM_178444</u> , <u>NM_178444.1</u> , <u>NM_178444.2</u> , <u>NM_178444.3</u> , <u>NM_178444.4</u> , <u>BB871615</u> , <u>BG292585</u> , <u>BY274501</u> , <u>NM_178444.5</u>
UniProt ID:	<u>Q9QXT5</u>
Synonyms:	VE-statin; Zneu1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)