

## Product datasheet for **MP204122**

### Ehd4 Mouse qPCR Primer Pair (NM\_133838)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	98878
Forward Sequence:	GCGAGAACTTATCTACAGGCTCC
Reverse Sequence:	GGTGAAGTCGTAGTTCTCCAGC
ACCN:	<u><a href="#">BC007480</a></u> , <u><a href="#">BC012272</a></u> , <u><a href="#">NM_133838</a></u> , <u><a href="#">NM_133838.1</a></u> , <u><a href="#">NM_133838.2</a></u> , <u><a href="#">NM_133838.3</a></u> , <u><a href="#">NM_133838.4</a></u>
UniProt ID:	<u><a href="#">Q9EQP2</a></u>
Synonyms:	2210022F10Rik; AI197390; AI846352; AV006278; Past2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)