

Product datasheet for **MP204140**

Elk1 Mouse qPCR Primer Pair (NM_007922)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13712
Forward Sequence:	GGAACAAGCTCTGGTCTTCAGG
Reverse Sequence:	GGACTCAGAGTGCTCCAGAAATG
ACCN:	BC054474 , NM_007922 , NM_007922.1 , NM_007922.2 , NM_007922.3 , NM_007922.4 , NM_007922.5 , BC005706 , BY728635
UniProt ID:	P41969
Synonyms:	Elk-; Elk-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)