

Product datasheet for **MP204185**

Exog Mouse qPCR Primer Pair (NM_172456)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	208194
Forward Sequence:	TTGTCGTGAGCTGACGGAGAGA
Reverse Sequence:	GCCACGTTGTCCTCACCAATCA
ACCN:	BC119024 , NM_172456 , NM_172456.1 , NM_172456.2 , NM_172456.3 , BC052171 , BC059920 , BC064730 , BC094589 , BC144729 , BI135003 , BY113415 , NM_172456.4
UniProt ID:	Q8C163
Synonyms:	AW557704; Endogl1; Endogl2; EngI; ENGL-a; ENGL-B; Engla; Englb
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)