

Product datasheet for **MP204262**

Epha2 Mouse qPCR Primer Pair (NM_010139)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13836
Forward Sequence:	GGCTGTACTCAAGTTTACCACCG
Reverse Sequence:	CCGCTTTCAGTGTCTTGATGGC
ACCN:	<u>NM_010139</u> , <u>NM_010139.1</u> , <u>NM_010139.2</u> , <u>NM_010139.3</u> , <u>BC140960</u> , <u>BB817901</u> , <u>BC006954</u> , <u>BY024029</u>
UniProt ID:	<u>Q03145</u>
Synonyms:	AW545284; Eck; Myk2; Sek-2; Sek2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)