

Product datasheet for **MP204318**

Ergic2 Mouse qPCR Primer Pair (NM_026168)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67456
Forward Sequence:	TCACTGACTCCAGATGCGTGCA
Reverse Sequence:	TGGTTGACAAGTGCTGCCAAATG
ACCN:	BC006895 , BC018188 , BC026558 , BC064749 , NM_026168 , NM_026168.1 , NM_026168.2 , NM_026168.3 , NM_026168.4
UniProt ID:	Q9CR89
Synonyms:	1200009B18Rik; 4930572C01Rik; AA408438; Ptx1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)