

Product datasheet for **MP204322**

Ero1b Mouse qPCR Primer Pair (NM_026184)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67475
Forward Sequence:	CTGGATGACTGTGAGCAGGCTA
Reverse Sequence:	TCACAAAAGTGGTCCTGCGAATC
ACCN:	<u>BC058721</u> , <u>BC111900</u> , <u>BC112426</u> , <u>NM_026184</u> , <u>NM_026184.1</u> , <u>NM_026184.2</u> , <u>BC145888</u> , <u>BC145890</u>
UniProt ID:	<u>Q8R2E9</u>
Synonyms:	1300013B24Rik; 1700065B09Rik; AI447560; ero1-beta; Ero1b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)