

Product datasheet for **MP204327**

Epb41I3 Mouse qPCR Primer Pair (NM_013813)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13823
Forward Sequence:	TTGAACGCTCGTCCAGCAAACG
Reverse Sequence:	AGAGATGCCCTTTGTTGTGGCG
ACCN:	<u>NM_013813</u> , <u>NM_013813.1</u> , <u>BC079858</u> , <u>BC152368</u> , <u>BC152548</u> , <u>BC167168</u> , <u>NM_013813.2</u>
UniProt ID:	<u>Q9WV92</u>
Synonyms:	4.1B; Dal1; DAL1P; Epb4.1I3; NBL3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)