

## Product datasheet for **MP204331**

### Dmtn Mouse qPCR Primer Pair (NM\_013514)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13829
Forward Sequence:	CATCCTGGACATTGAGCGACCT
Reverse Sequence:	GGGATGTGGATTTGGGTGACAG
ACCN:	<u><a href="#">BC016897</a></u> , <u><a href="#">BC037021</a></u> , <u><a href="#">NM_013514</a></u> , <u><a href="#">NM_013514.1</a></u> , <u><a href="#">NM_013514.2</a></u> , <u><a href="#">NM_013514.3</a></u> , <u><a href="#">NM_013514.4</a></u> , <u><a href="#">BY000876</a></u> , <u><a href="#">BY270342</a></u> , <u><a href="#">BY282941</a></u> , <u><a href="#">BY601680</a></u> , <u><a href="#">BY789340</a></u>
UniProt ID:	<u><a href="#">Q9WV69</a></u>
Synonyms:	AI325486; dematin; Epb4.9; Epb49
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)