

Product datasheet for **MP204350**

Ebag9 Mouse qPCR Primer Pair (NM_019480)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55960
Forward Sequence:	AAGATGCAGCCTGGCAAGCAGA
Reverse Sequence:	TCTGCTCCTTCTTCATCAGGCG
ACCN:	<u>BC108397</u> , <u>BC119118</u> , <u>BC119120</u> , <u>NM_019480</u> , <u>NM_019480.1</u> , <u>NM_019480.2</u> , <u>NM_019480.3</u> , <u>NM_019480.4</u> , <u>BC026717</u> , <u>BC144740</u> , <u>NM_019480.5</u>
UniProt ID:	<u>Q9D0V7</u>
Synonyms:	AI835379; Rcas1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)