

## Product datasheet for **MP204356**

### Ei24 Mouse qPCR Primer Pair (NM\_007915)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13663
Forward Sequence:	CTTGGCTTTCCTCACAGCAATGC
Reverse Sequence:	GGAGTCTTTGCTTCATTGGCGC
ACCN:	<a href="#">NM_007915</a> , <a href="#">NM_007915.1</a> , <a href="#">NM_007915.2</a> , <a href="#">NM_007915.3</a> , <a href="#">NM_007915.4</a> , <a href="#">NM_007915.5</a> , <a href="#">BC140389</a> , <a href="#">BC146473</a> , <a href="#">BU148129</a> , <a href="#">BY282890</a>
UniProt ID:	<a href="#">Q61070</a>
Synonyms:	AA536736; AI115355; PIG8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)